Science. Applied to Life

Patients Vary.Our Standards Don't.

Healthcare Associated Infections (HAI's) are the most common complication affecting patients in Hospitals.¹

Patient safety and best practice start in the CSSD.

There are five essential steps in sterile processing:

1) Clean, 2) Prep & Pack, 3) Sterilise, 4) Store, and 5) Use.

The threat of SSIs caused by non-sterile devices makes following every step in the sterilisation process absolutely critical.



Healthcare Associated Infections (HAI's) are the most common complication affecting patients in Hospitals.

Each Year around 200,000 HAI's are contracted by patients in Australia.

Reference: NSQHS Standard.

How do we manage this risk? The Australian Commission on Safety and Quality in Healthcare NSQHS standards state that "Reprocessing of equipment and instrumentation meets current best practice guidelines."

Total national number of bed days due to **Surgical Site Infections** for a one year period was estimated to be 206,527 bed days.

Reference in the NSQHS Standard 3

NSQHS 3.16.1 states "Compliance with relevant national or international standards and manufacturer's instructions for cleaning, disinfection and sterilisation of reusable instruments and devices is regularly monitored"



What does the international community say about managing the Risk of SSI in relation to CSSD?

Using National and International standards as suggested in the Australian Commission on Safety and Quality in Healthcare, what does Best Practice look like?

Criteria for Load Release	AS/NZS 4187:2014	ANSI/ AAMI ST79: 2013	AORN 2016	APSIC 2017	CDC 2008	World Health Organisation (WHO) 2016	Best Practice
Equipment							
Temperature	Every Load	Every Load	Every Load	Every Load	Every Load	Every Load	Every Load
Pressure	Every Load	Every Load	Every Load	Every Load	Every Load	Every Load	Every Load
Time	Every Load	Every Load	Every Load	Every Load	Every Load	Every Load	Every Load
Bowie & Dick test	Daily	Daily	Daily	Daily	Daily	Daily	Daily
Load							
Biological Indicator (BI) in test pack	Optional	Weekly, Preferably Daily	Weekly, Preferably Daily	Minimum Daily	Weekly, Preferably Daily	Minimum Daily	Weekly, Preferably Daily
Biological Indicators - Implants test pack		Every Load	Every Load	Every Load	Every Load	Every Load	Every Load
Chemical Indicator (CI) in test pack	Optional	Optional	Optional			Every Load	Every Load
Pack							
Internal Chemical Indicators	Optional	Every Pack	Every Pack	Every Pack	Every Pack	Every Pack	Every Pack
Exposure							
External Chemical Indicators	Every Pack	Every Pack	Every Pack	Every Pack	Every Pack	Every Pack	Every Pack

"A Quality Monitoring Program that includes **Physical monitors** (eg printouts, digital readings, graphs, gauges), Chemical Indicators, and Biological **Indicators** should be used to verify that conditions necessary for steam sterilisation have been met".

Reference: VI.c. (p827) AORN 2016 guidelines (Strong evidence).

Monitoring Critical Variables in the Ste

Equipm



Who gets the results?

CSSE

Possible reasons why indicators may fail in load or pack monitoring

Incorrect pack orientation

Too many packs in the load

Placing packages too close to each other

Incorrect cycle selection

Inadequate steam quality

Sterilisation process monitoring devices include Physical monitors, **Chemical Indicators and Biological Indicators.**

> Each of these devices plays a distinct and specific role in sterilisation process monitoring, and each is indispensable to sterility assurance [AAMI ST79:2013]

erilisation process						
ent	Load	Pack				
ce in liser cy	Indicates pack contents exposed to correct conditions	Confidence that pack contents exposed to correct conditions				
	\bigotimes	\bigotimes				
	\bigotimes	\bigotimes				
	\bigotimes	\bigotimes				
	\checkmark	N/A				
ick	In test pack	\checkmark				

)	CSSD	Theatre

A reminder on why this is important. What if you do get a BI failure?

A failure in one of the measures of the sterilisation process would require your Recall Procedure to be put into place immediately, to recall every load since the last negative BI result.*















How is your patient safety compromised?

*Based on a CSSD that processes 5 loads per day





What does the **Risk Profile** in your CSSD look like today?

Criteria for Load Release	Best Practice	Current Practice today	What do you want to do?	By when?			
Equipment							
Temperature	Every Load						
Pressure	Every Load						
Time	Every Load						
Bowie & Dick test	Daily						
Load							
Biological Indicator (BI) in test pack	Weekly, Preferably Daily						
Biological Indicators - Implants test pack	Every Load						
Chemical Indicator (CI) in test pack	Every Load						

Pack

Internal Chemical Indicators	Every Pack			
------------------------------------	------------	--	--	--

Exposure

External Chemical Indicators	Every Pack			
------------------------------------	------------	--	--	--



Today, someone's immediate family member is visiting the operating theatre. Patient Safety in Surgery starts in the CSSD.

How can 3M help?

Control	Frequency	Reference ISO standard	Available 3M Product		
Equipment	Daily	Daily ISO 11140-4 ISO 11140-3	Comply 1300 Bowie Dick Test Pack		
			Comply 1301 Bowie Dick Test Pack		
			4208 ETS		
Load	Daily	ISO 11138-1 ISO 11138-3	Attest™ 490 Super-Rapid Auto-Reader (1 hour result)		
			1492V Biological Indicaators (box 50)		
				41482V Biological Indicator + Type 5 Sterigage Integrator Test Pack (24 Test Packs + 24 Controls)	
			Attest [™] 390 Rapid Auto-Reader		
			1291/1292 Biological Indicators (box 50) 1296 Rapid Biological Test Pack (25 test packs + 25 controls)		
Pack	Every Pack	ISO 11140-1	Type 5 Comply Sterigage Steam Chemical Integrators (1243A)	120 Court - Sendlar - Courts	
Exposure	Every Pack	ISO 11140-1	Comply Tape 1322-12mm	(A)	
			Comply Tape 1322-18mm	er er	
			Comply Tape 1322-24mm	4	
			Comply Tape 1355-18mm		
			Comply Tape 1355-24mm	AND -	

Always read the label and use only as directed



3M Infection Prevention Division

3M Australia Pty Limited Bldg A, 1 Rivett Road North Ryde NSW 2113 1300 363 878 www.3M.com.au/healthcare

References:

- The Australian Commission on Safety and Quality in Healthcare NSQHS Standard 3 2012: factsheet
 AS/NZS 4187:2014 Reprocessing of reusable medical devices in health service organisations
 Centers for Disease Control and Prevention, Guideline for Disinfection and Sterilisation in Healthcare Facilities, 2008
 The APSIC guidelines for disinfection and sterilisation of instruments in health care facilities, 2010
 S. World Health Organisation (WHO), Decontamination and Reprocessing of Medical Devices for Health-care Facilities 2016
 ANSI/AAMI ST79:2013, Comprehensive guide to steam sterilization and sterility assurance in health care facilities
 7. Association of periOperative Registered Nurses (AORN), Guidelines for Perioperative Practice, 2016 Edition.